



New Optical Thermal Analysis Methods: Absolute, Differential, TGA, Atmosphere

Expert Lab Service is a boutique of ceramic engineering solutions, materials analysis services and **tailor-made laboratory instruments.**

Outline

- Company presentation
- The new **Absolute Optical Platform**
- The new **1700°C MoSi₂** furnace for **Vertical Dilatometer + Microscope**
- The new **Differential** Dilatometry and Fleximetry
- Simultaneous **ThermoGravimetry**
- **Controlled Atmosphere**: vacuum and inert
- **Modelling** software

What we do

- **Measurement instruments** for optical thermal analysis: ELS-MDF (microscope, dilatometer, fleximeter)
- **Ceramics Genome:** Ceramic laboratory information management and modeling software
- **Consulting and training**
- **Laboratory tests**

Selected customers

→ CERAMIC TILES MANUFACTURERS

LAMINAM



→ RAW MATERIALS SUPPLIERS



→ RESEARCH INSTITUTIONS



→ GLAZE PRODUCERS



TORRECID



esmalglass.itaca
grupo

COLOROBBIA
ITALIA

→ TECHNOLOGY SUPPLIERS



→ ADVANCED CERAMICS, AEROSPACE



+GF+ PRECICAST



Optical Thermal Analysis



Optical thermal analysis instruments are able to measure thermophysical properties of materials **without any contact between the measurement system and the sample.**

They all work by applying **computer vision techniques** to the sample's image (or some portions of it).

Measurable properties: coefficient of thermal expansion, melting point, elastic tensions, plastic deformation, contact angle, surface tension, viscosity estimation.


Seeing > Touching

Being contactless, we avoid most of the interference of the instrument on the sample.



Kinetics > Equilibrium

The instrument does not heat up during the measurement. This allows to observe the sample behavior with very fast heating rates.



Materials > Standards

Standards can be too vague to be useful, or of little applicability for your material. We always prioritize materials over standards.

From ISO 540 on ash fusibility, 1995:

3.1 deformation temperature (abbreviation DT):

The temperature at which the *first signs of rounding*, due to melting, of the tip or edges of the test piece occur.


From ASTM E381-19 on CTE with thermomechanical analysis:

6.1.3 Sensing element... $\pm 50\text{nm}$ resulting from changes in length of the specimen.


6.3 Micrometer... with a range up to 10mm to determine specimen dimensions within **$\pm 25000\text{nm}$** .

Optical Thermal Analysis for Ceramics

- We measure the deformation that heat and temperature causes in a material, across these behaviours:






Design of an optical thermal analysis instrument




- The optical system brings the *image* of the sample on camera.
- Resolution is limited by the optical system: max ~0.5 microns.
- We can focus a part of the sample and follow its **displacement**, or its entirety to follow its **morphing**.

Optical Measurements: Displacements










- **Displacement:** maximum zoom ($1\text{px}=0.4\mu\text{m}$) on the smallest possible segment of sample border. We get a (sometimes not so) straight line.
- Position through linear regression + mean line-line distance + statistical filtering
- Delta with initial position = displacement
- **Single point: Fleximeter**
- **Two points: Horizontal/Vertical Dilatometer**
- **Three points: Absolute Fleximeter**




Optical Measurements: Heating Microscope

Follows the change of the entire silhouette of the sample, as it melts.




Microscope, Dilatometer, Eleximeter




ELS-MDF

The all-in-one optical thermal analyzer




State of the Art: Motorized Optics



- Full HD cameras (1280x920 px) magnified to 0.5 μ m/px are able to follow 640 μ m, or **1.3%** of expansion/sintering over 50mm.
- When the sample is going to exit this window, a motor moves the camera, takes a new zero and restarts the measurement.
- Drawbacks:
 - Each movement causes 0.5 μ m of error.
 - The dynamics occurring during the intervention time is lost
 - The initial length must be measured with an external micrometer, usually having 20 times less resolution than the instrument itself (10 μ m against 0.5 μ m!)

Next Iteration: Absolute Optical Platform




High
Precision
Optics (HPO)

High Range
Optics (HRO)

- The **Absolute Optical Platform** features a total of **60Mpx** of sensors and a massively parallel computer vision software to simultaneously decompose and analyze multiple huge frames.
- The entire measurement is executed **without any camera movement**.
- The position of the optics relative to each other is fixed in a robust **Invar** assembly, **immune to thermal expansion**.
- The field of view is big enough to make possible, for the first time in a commercial instrument, **Differential Optical Measurements**.
- **Initial dimension measurement** is possible.
- **No complex alignments. Easier maintenance. Self-installation.**
- The entire optical bench can translate XY in order to execute **8 functions**


Absolute Optical Platform: High Range



High
Range
Optics

- Two **high range optics** (HRO) for:
 - HRO **Horizontal Dilatometer**, to follow **sintering** up to **35%** with 3.4um resolution over 50mm (70ppm).
 - HRO **Fleximeter**, for **pyroplasticity**, up to 6mm.
 - **Heating Microscope** up to 3 samples (variants up to 4 and 8), up to 10x10mm
 - HRO **Differential Horizontal Dilatometer**, to load non-sticky samples (see later).

Absolute Optical Platform: High Precision




**High
Precision
optics**

- Two **high precision optics** (HPO) for:
 - HPO **Dilatometer** to follow **CTE** up to **4%** with 0.5um of resolution over 50mm (10ppm).
 - HPO **Absolute Fleximeter**: One HP optic in the middle and two HR at the sides, to execute the circular approximation and the Absolute Fleximeter function.
 - HPO **Differential Fleximeter**: by framing both a non-bending reference and the sample.
 - **HP Microscope** samples of less than 1x1mm.

→

Absolute Optical Platform: 9 functions




High Precision Optical (HPO):

- 1) Horizontal Dilatometer (CTE, 10ppm res)
- 2) Differential Horizontal Dilatometer (non-sticky)
- 3) Differential Fleximeter (<150um tensions)
- 4) Absolute Fleximeter (>150um tensions)
- 5) Microscope (max 1.5mm height)

High Range Optical (HRO):

- 6) Dilatometer (for sintering, 70ppm res)
- 7) Differential Dilatometer (or double sample)
- 8) Fleximeter (for pyroplasticity)
- 9) Microscope (x3, x4, x8, max 10mm height)


Absolute Optical Platform: Vertical 1700°C



High
Precision
Optics (HPO)

High Range
Optics (HRO)

- MoSi2 are less fragile, less expensive and easier to replace than Pt/Rh elements
- **Heating Microscope:**
 - The vertical design of the furnace allow to run sintering tests while framing the entire sample (15mm of height)
 - A dedicated, microscope-only design can host x8 samples
- **Differential Vertical Dilatometer:**
 - 15mm sample, 0.5um resolution, for accurate CTE




MoSi2 Furnace


The Optical Differential Concept

- Even optical instruments cannot be immune to **thermal anisotropies** affecting the sample, the furnace and the thermocouple. Fast heating rates, possible because of the cold measurement system, exacerbate this issue - especially at lower temperature (<400°C).
- In the **vertical** dilatometer configuration, moreover, temperature gradients in air moving parallel to the measurement direction can cause **optical illusions and aberrations**.
- Thanks to the new high resolution sensor, **both the sample and a calibration reference can be framed** simultaneously.
- 3 variants:
 - **Differential Fleximeter**: a non-bending reference and the sample are stacked one above the other, with a tiny space between them where the sample can bend. **One camera** frames the reference and the sample at the same time.
 - **Differential Vertical Dilatometer**: sample and reference stand close to each other. Their tips, separated by a tiny space, are framed by the **same camera**.
 - **Differential Horizontal Dilatometer**: sample sits above the reference (which can be an alumina plate). Two cameras frame the interface point between sample and reference.


The Optical Differential Fleximeter




The Optical Differential Vertical Dilatometer



- Traditionally, Optical Vertical Dilatometer reads a bad signal at low T, showing an **apparent shrinking**. Caused by difference in air motion between the top and the bottom of the sample.
- The situation was better with small furnaces (impractical with MoSi₂). Was used only for **sintering** studies.
- The **Differential** technique **measures the illusion on the reference** and subtracts it from the sample.
- Leading to **good CTE** measurements in Vertical mode even at low T.



The case for Absolute+Differential Fleximeter



- Measuring in three-points gives an additional information about the material: is it bending like an **elastic** or a **plastic** material?
- Elastic bending follows **Timoshenko beam theory**. The material curves along a **circular** profile with a curvature radius depending on elastic moduli and CTE mismatch. Even a uniform a material will bend elastically if its modulus drops with rising T.
- **Pyroplastic** bending follows the parabolic profile of bending moments.
- Absolute+Differential can measure:
 - The relative modulus change in the elastic phase
 - The exact onset of the transition elastic→plastic

$$r = \frac{L}{6} \frac{3(1+m)^2 + (1+mn)(m^2 + \frac{1}{mn})}{(\alpha_2 - \alpha_1) \Delta T (1+m)^2}$$

m = thickness ratio, L_1/L_2

L = total thickness, $L_1 + L_2$


n = moduli ratio, E_1/E_2

α = coefficients of thermal expansion

ΔT = temperature delta, causing expansion of each layer

r = curvature radius, produced by the differing expansion


Optical Differential Horizontal Dilatometer



- The **Differential Horizontal Dilatometer** techniques allows to cancel the temperature effect. The reference sample has a large contact surface with the sample, lower thermal mass (thin plate) and can also have higher thermal diffusivity (in case of metal plate).
- This function can also measure two independent samples without any differential calculation. Doubles the productivity.
- Caveat: sample must be not stick with reference or with the second sample during the measurement!



Simultaneous Thermogravimetry: Dil-TGA, Flex-TGA

- Combined techniques offer the obvious advantage of **saving time**, as two or more properties are measured during the same test run.
- The most interesting point is **commensurability**:
 - Typical TGA setup analyzes **few mg** of material in a **tiny and slow furnace**
 - While the results can be extremely accurate, they might be not much **representative of what happens** in a **bigger** sample at **faster** rates (eg: rate of gas exchange)
- As dilatometry was traditionally a contact measurement, **no simultaneous TGA was conceivable!**
- Simultaneous TGA allows to measure weight change on the **same sample**, in the **same furnace**, with the **same thermal cycle** and at the **same time** as the dilatometric or fleximetry signal.
Commensurable!
- We mounted the entire sample holder assembly on different types of load cells (here *one* of the iterations).





Thermogravimetry

- We are testing load cells in the range 1-10mg. No ambition to compete with dedicated TGA: we aim at industrial ceramics, weight losses in the order of 2-5%. Aim at 0.1% resolution.
- Even with 10mg resolution, overall results over a typical 3gr sample are comparable with TA Instrument's QSeries DSC-TGA.
- Not aiming at measuring during Microscope analysis. Samples are too small.
- Delivery expected in Q2 2024.



Controlled Atmosphere and Vacuum

- The entire furnace, sample holders, heating elements contacts and thermocouple is enclosed in a vacuum chamber. Mild vacuum down to 5mBar.
- Mounting our standard platinum heating elements. Tested up to 1300°C in vacuum on spacecraft thermal shields.
- By purging with inert gas and then filling slightly above atmospheric pressure, ensures absence of O₂.
- Tested over 3D additive manufacturing metal alloys without signs of oxidation
- Delivery expected in Q4 2024.



Ceramics Genome: *bring ceramics alive!*

- Ceramic modeling software: from the composition of raw materials and minerals to the laboratory measurements.
- A **vertical solution specialized** in ceramics laboratory information management. Entirely web based.
- Launched for beta testing in 2022, we are now customizing the user experience.
- Two main components: Composition and **Measurements**.
- The Measurements components creates web Reports easy to **search and view**.
- Our target is to provide **standardized** plotting and calculation tools covering **most ceramists' needs** in a **simplified** package.

Date	O...	Name	T...	Duration	°C	°C/min
9/20/2023, 8:15:39 AM		Au granitico 20230920		78.84	1082.49	50.00
9/19/2023, 4:31:00 PM		ag granitico		25.54	979.59	50.00
6/30/2023, 11:59:52 AM		oro platino 1450		28.44	1401.48	50.00
6/30/2023, 11:59:52 AM		Alloy Au Pd 1450		28.44	1401.48	50.00
6/30/2023, 10:34:14 AM		Alloy Au Pd		53.57	1401.39	50.00
6/30/2023, 8:07:25 AM		Alloy Au Pd		142.56	742.51	5.00
6/29/2023, 1:54:50 PM		Alloy Au Pd		72.43	1398.74	50.00
5/23/2023, 7:46:48 AM		Indium 20230523 -20 80 - 1 200		83.73	160.86	20.00
5/22/2023, 8:15:25 AM		Indium 20230522 -20 80 - 1 200		83.70	160.69	20.00
5/19/2023, 7:49:36 AM		Indium 20230519 -20 80 - 1 200		86.20	163.26	20.00
5/18/2023, 1:39:29 PM		Indium 20230518 -20 80 - 1 200 ii		97.55	174.67	20.00
5/18/2023, 1:39:29 PM		Indium 20230518 -20 80 - 1 200 ii		97.55	174.67	20.00
5/18/2023, 10:45:03 AM		Indium 20230518 -20 80 - 1 200		104.87	181.90	20.00
5/18/2023, 10:45:03 AM		Indium 20230518 -20 80 - 1 200		104.87	181.90	20.00
5/18/2023, 8:01:15 AM		Indium 20230518		38.82	216.49	5.00
5/17/2023, 2:15:28 PM		Saccarin dried 20230517 stasi iii		40.77	249.50	20.00
5/17/2023, 11:12:23 AM		Saccarin dried 20230517 stasi ii		40.77	249.55	20.00
5/17/2023, 8:28:04 AM		Saccarin dried 20230517 stasi		40.84	249.38	20.00
5/16/2023, 4:27:10 PM		Saccarin dried 20230516 stasi ii		40.67	249.54	20.00
5/16/2023, 2:49:58 PM		Saccarin dried 20230516 stasi		40.60	249.49	20.00
60		60		65.26	26.63	


Download table Load all

< delta t vetro finestra, HSM

📄 🔍 🔄 ⚙️

Main Logging File Bill Sample0 Sample1 Sample2 Admin

Name basso sx
Sintering 693.67 °C (3956 s)
Softening 804.36 °C (4621 s)
Sphere 896.42 °C (5173 s)
Half Sphere --- °C (--- s)
Melting --- °C (--- s)




I: Sintering
R: Softening
S: Sphere

Shapes Images

Show


Ceramics Genome: Overlay curves

- Select any number of tests to create an overlay document that can be saved.
- Each curve can show characteristic points if found during the test.




Ceramics Genome: Coupling temperature

- Detect the coupling temperature from body and glaze dilatometric curves and the product fleximeter curve.
- Models the theoretical Timoshenko curvature from the dilatometer tests and detects the point at which the real fleximeter departs from the prediction as the coupling T.



Ceramics Genome: Viscosity

- Based on the equations from Adcock, Drummond, McDowall (1959)
- Measurement Inputs: Flex, to calculate the speed of deformation; Dilatometer, to correct the volume of the sample; TGA, to correct the density. Can use Flex-TGA and Dil-TGA signals.
- Needs initial density and sample's geometry if not saved by the operator in the test file.




Ceramics Genome: Quartz estimation

- Estimate the volumetric quartz content based on the shape of the quartz transition.
- Select the curve, input bulk density if known, play.
- Outputs the estimated volume and weight, and the model plot showing the discrepancy between the predicted curve and the real one.



Main Plot Data Output Admin

Success	×
Message	---
Quartz Volume	11.09 %
Quartz Weight	12.56 %
Finest quartz fraction	---
Intermediate quartz fraction	---
Rough quartz fracion	---





Thank you

Expert Lab Service is a boutique of ceramic engineering solutions, materials analysis services and **tailor-made laboratory instruments.**